



ENH 286 -

Merial

Honeywell 99EXhc Mobile
Computer Training

Author:	Letitia Sims
Creation Date:	June 23, 2014
Last Updated:	June 23, 2014
Document Ref:	Merial_5174
Version:	1.0




Table of Contents

A.	Keyboard	4
B.	Battery	8
C.	Sign into Scanner (Mobile Computer)	10
D.	Creating Labels:	12
E.	Planting Process	18
F.	Harvest Process	24
G.	Exception process for vHVT virus	36
H.	Maintenance Reason Codes:	44
I.	Edit Control Report	54
J.	Requesting AI Control Report from Oracle Form:	59
K.	End of Process	63












A. Keyboard

55-Key Full Alpha/Numeric Keyboard












Note: The flashlight function () is disabled on all 99GX models.

Using the Function Keys

Name	Key	Function
Backlight		Turns the keyboard backlight on and off.
Flashlight	 & 	Turns the flashlight on and off. The flashlight is located on the back panel of the terminal. <i>Note: The flashlight function is disabled on all 99GX models.</i>
Backspace (BKSP)	 55-Key	Moves the cursor back one space each time the key is pressed. If you are typing text, it deletes the previous character each time it is pressed.
Delete (DEL)	 55-Key	Deletes the next character each time the key is pressed. This key only appears on the 34-key and 55-key keypad options.
Enter (ENT)		Confirms data entry.
Escape (ESC)		Cancels the current action.
Power Key		Puts the terminal in and wakes the terminal from Suspend Mode; see Suspend Mode on page 2-11.
SCAN Key		Activates the scan and "wakes" the terminal from Sleep Mode. Its position allows convenient one-handed image-taking and/or bar code decoding.
Space (SP)	 55-Key	Moves the cursor one space.
Tab		Moves the cursor to the next tab stop or the next control (on a form).

Using the Navigation Keys

Located in the center of each keyboard for easy access with either hand, the navigation keys navigate the cursor through application screens.

Key	Function
	Moves the cursor up one row or line.
	Moves the cursor down one row or line.
 & 	Raises the volume.
 & 	Lowers the volume.
 & 	Moves the cursor up one page.
Key	Function
 & 	Moves the cursor down one page.
	Moves the cursor one character to the right.
	Moves the cursor one character to the left.

Note: Additional functionality varies according to the application in use.

Additional function keys:

 +  exit current screen

 +  return to previous screen

 +  list of values

 +  soft reset (warm boot)

Power key   :

- turns on device
- puts the device in “suspend” mode (i.e. sleep)

B. Battery


Removing the Main Battery Pack

When removing a battery from the terminal, put the device in **Suspend Mode** (see page 9) before removing the battery door. Once the battery door is removed, **wait at least 3 seconds** before removing the main battery. This process allows the device to shut down properly and to maintain memory during the battery swap.

Connect the terminal to one of the 99EX series charging peripherals to charge the main battery.

The charging time for the main battery pack is 4 hours for the standard battery or 6 hours for the extended battery.

Honeywell recommends charging the Dolphin terminal for at least 24 hours prior to initial use to ensure the internal backup battery is fully charged.





Timeout settings:


IMPORTANT --

Each mobile device is set to 60 minutes of connectivity to the network. So long as the user continues to use the device under 60 minutes, the data entered will not be lost. All data 'submitted' will be captured. The device will go into a "sleep" mode, meaning the screen will go black, but this does not mean that the connectivity is lost.

After 60 minutes, the user will have to sign back into the device with the Oracle userid and password; however, any data that was not 'submitted' will be lost and the user will have to re-enter the data.

C. Sign into Scanner (Mobile Computer)


- User clicks  Power key
- (Note: if the Naurtech screen returns, click on **"Session"** → **"Connect"**)
- User clicks  or chooses '1' for Default




- Enter your Oracle user name and password (you must have an Oracle userid and password in order to access this program)



- Password will be show on the screen as asterisks



- The main menu will appear:




D. Creating Labels:


Creating Rack Code Labels:

- Users will be provided with a Microsoft Excel Marcos label utility that will provide the ability to create Rack Code labels.
- Enter the Rack Code and the number of that label desired

The screenshot shows an Excel spreadsheet with a grid from column B to N and row 5 to 36. A form is overlaid on the grid, containing the following elements:

- A text box labeled "Rack Code" with the value "1234IJKL" entered.
- A text box labeled "Number of Lables to Print" with the value "1" entered.
- A small empty text box below the "Number of Lables to Print" field.
- Two buttons at the bottom: "Clear Fields" and "Label Print".

- Press 
- The screen will automatically go to "Sheet 2" and ask the user if they want to "save the changes made on the label.csv?"
- The user must choose either 'Yes' or 'No'. If you choose cancel, you will receive a debug error message and your labels will not print.



- The Rack Code labels will be printed from the Zebra S4M label printer.




- These are permanent labels used to mark the racks once, unless otherwise repeated.


Creating Ser Code (Lot) Labels:

When starting a process for a cell start, the user must first create Ser Code (lot) barcode labels.

- Choose Request lot labels and the Print Unique Labels screen will return




- User will choose the org (organization) code
- Enter **GMS** (if requested)




Each rack used in the process of this cell start will be associated with the Plant Number and Ser Code (lot) label based on the number of racks planted.


- User will know the cell start (lot number) and the number of racks they will need for the process.
- (**Note:** Lot number must be in capital letters. Use soft keyboard with stylus and set it to CAPS. Then close the soft keyboard).




- The racks will have permanent barcode labels on them and will be known as the "Rack Code" label


- User will manually enter the cell start Lot Number, Plant Number and number of labels (which will equal to the approximate number of racks to be planted)




- User presses the  button to highlight "Submit" and screen returns to main menu
- If user enters wrong information, user can arrow down to choose to "Cancel" and exit the screen to start over.





- User has to enter a lot number that is for a Cell Start item or the error message should be "No Result is found"




- Once submitted, the Ser Code (lot) labels are printed on the Zebra S4M label printer
- (**Note:** the barcode is a unique identifier that represents the lot number and the plant number to allow both to be captured with one scan):

Lot# CEF12141C		Plant# PLANT5	
			
Tank#	Rack#	Hood#	Tech
T—	R—	H—	

Lot# CEF12141C		Plant# PLANT5	
			
Tank#	Rack#	Hood#	Tech
T—	R—	H—	

Lot# CEF12141C		Plant# PLANT5	
			
Tank#	Rack#	Hood#	Tech
T—	R—	H—	

Lot# CEF12141C		Plant# PLANT5	
			
Tank#	Rack#	Hood#	Tech
T—	R—	H—	

When there is more than one item number with the same cell start lot number:

If there is a lot number that is used by two item numbers, you will be prompted to choose which item you want:


```

C:\ Telnet 10.39.0.206
Print Unique Labels
Lot Numbe>CEF01031A
Plant Num:
# of Labe:
<Submit>
<Cancel>


```

- The following screen will return


```
C:\ Telnet 10.39.0.206
Lot Number!Item Co
1<CEF01031A !2455000
2<CEF01031A !245500>
```

- Use your right arrow cursor  to scroll to the right in order to see the full item number

```
C:\ Telnet 10.39.0.206
ot Number!Item Code
1<01031A !24550000 >
2<01031A !2455000B >
```

- Choose the desired item with the up or down cursor or the number choice and press  choice and press


E. Planting Process

Bottles are filled and placed on the roller racks. This next section is how the scanner will be used during the vaccine planting process.


Process Control → Process AI (screen)

The process control screen associates the cell start lot number, plant number, virus, and process to an available rack(s) by entering the rack and scanning the Ser label (lot) barcode's unique number. This is done during the planting process. It also establishes the control number for reconciling remaining processes to the number of racks planted.


- User chooses Process control screen from main menu




- User manually enters the Lot Number and Plant Number




- User clicks  and chooses the Lab from the list of values



- User clicks  and chooses the Virus from the list of values



- If HVT, Rispsens, SB-1 Seed or SB-1 Serial is chosen, the system will default to "Plant & Infect" for the first process



- If vHVT (Vaxxitek) is chosen, the system will default to "Plant" for the first process


```

C:\ Telnet 10.39.0.206

Process AI
Lot Numbe>CEF12141C
Plant Num:PLANT5
Lab      >TC1
Virus    >HUT <Uaxx
Process  >1-Plant
Rack Code:
Ser Code :
<Continue>
<Submit>
<Cancel>

```

- User clicks  to scan the Rack Code Label




- Remaining Rack Codes used in this example are:
- 1234EFGH
- 1234IJKL
- 1234MNOP

```

C:\ Telnet 10.39.0.206

Process AI
Lot Numbe>CEF12141C
Plant Num:PLANT5
Lab      >TC1
Virus    >HUT
Process  >1-Plant &
Rack Code:1234ABCD
Ser Code :
<Continue>
<Submit>
<Cancel>

```


- User click  to scan the Ser Code (lot) label (**note:** the unique number under the barcode will appear on the screen after the scan):

```

C:\ Telnet 10.39.0.206
Process AI
Lot Numbe>CEF12141C
Plant Num:PLANT5
Lab      >TC1
Virus    >HUT
Process  >1-Plant &
Rack Code:1234ABCD
Ser Code :25_
<Continue>
<Submit>
<Cancel>

```

- Label scanned was:

Lot# CEF12141C		Plant# PLANT5	
			
Tank#	25	Rack#	Hood# Tech
T —		R —	H —

Note: If user scans a Rack Code that is already associated with a different Ser Code (lot) label (other than what is currently being used), the error message should return as “Label for another lot”. Refer to business process rules and appropriate notifications.


```

C:\ Telnet 10.39.0.206
Process AI
Lot Numbe>CEF12141C
Plant Num:PLANT5
Lab      >TC1
Virus    >HUT
Process  >1-Plant &
Rack Code:1234ABCD
Ser Code :2
<Continue>
<Submit>
<Cancel>


Label for other lot


```

- When ready, the user can continue to associate the racks used with this lot:



- **Note:** If a Ser Code (lot) label has already been associated with another rack, the error message returns as "Rack already assoc"



- When all desired racks have been entered for association with the Ser Code and Rack Code, the user can arrow down  to choose "Submit"

```
C:\ Telnet 10.39.0.206
Process A1
Lot Numbe>CEF12141C
Plant Num:PLANT5
Lab      >IC1
Virus    >HUT
Process  >1-Plant &
Rack Code:
Ser Code :
<Continue>
<Submit>
<Cancel>
```

- The screen will return to the main menu

F. Harvest Process


The same Process AI screen will be used for the all of the remaining processes for the virus chosen.

- For the HVT, Rispens, SB-1 Seed and SB-1 Serial viruses, the only other process that is available after “Plant & Infect” is “Harvest”. However, “Harvest” will not be available in the list of processes if the “Plant & Infect” process has not been completed first.
- For the vHVT (Vaxxitek) virus, the first process available will be “Plant”. Once the “Plant” process has been completed, the next process will be available in the process list. This is “Infect Seed”, followed by “Harvest Seed”, “Infect Serial” and “Harvest Serial”. If the proceeding process has not been completed, the user will not be able to choose the next expected process for this virus.


1. HVT, Rispens, SB-1 Seed and SB-1 Serial process:

Below represents the example of the Harvest process for the HVT virus which is the same process for Rispens, SB-1 Seed and SB-1 Serial viruses:

- The user will manually enter the lot number and plant number
- The user will choose the lab and virus



- The user would choose the “Harvest” process from the list of processes:



- The user scans the desired Rack Code and Ser Code

```

C:\ Telnet 10.39.0.206
Process AI
Lot Numbe>CEF12141C
Plant Num:PLANT5
Lab      >TC1
Virus    >HUT
Process  >2-Harvest
Performed:
Rack Code:1234ABCD
Ser Code :25_
<Continue>
<Submit>
<Cancel>

```

- The user presses  to highlight Continue

```

C:\ Telnet 10.39.0.206
Process AI
Lot Numbe>CEF12141C
Plant Num:PLANT5
Lab      >TC1
Virus    >HUT
Process  >2-Harvest
Performed:
Rack Code:1234ABCD
Ser Code :25_
<Continue>
<Submit>
<Cancel>

```

- The screen will return with the number of racks harvested versus the total number of racks eligible for the Harvest process for this particular Cell Start Lot and Plant numbers

```

C:\ Telnet 10.39.0.206
Process AI
Lot Numbe>CEF12141C
Plant Num:PLANT5
Lab      >TC1
Virus    >HUT
Process  >2-Harvest
Performed:1 / 3
Rack Code:
Ser Code :
<Continue>
<Submit>
<Cancel>

```

- The user will enter the next Rack Code and Ser Code for the Harvest process

- The user will continue to enter all Rack Codes and Ser Codes until all desired racks have been entered for the Harvest process

```

C:\ Telnet 10.39.0.206

Process AI
Lot Numbe>CEF12141C
Plant Num:PLANT5
Lab      >TC1
Virus    >HUT
Process  >2-Harvest
Performed:1 / 3
Rack Code:1234EFGH
Ser Code :26
<Continue>
<Submit>
<Cancel>

```

```

C:\ Telnet 10.39.0.206

Process AI
Lot Numbe>CEF12141C
Plant Num:PLANT5
Lab      >TC1
Virus    >HUT
Process  >2-Harvest
Performed:2 / 3
Rack Code:
Ser Code :
<Continue>
<Submit>
<Cancel>

```


```

C:\ Telnet 10.39.0.206


Process AI
Lot Numbe>CEF12141C
Plant Num:PLANT5
Lab      >TC1
Virus    >HUT
Process  >2-Harvest
Performed:2 / 3
Rack Code:1234IJKL
Ser Code :27
<Continue>
<Submit>
<Cancel>

```

- Once all racks have been entered for the Harvest process, the "Performed" count will be equaled.



- Use Page Down button  and highlight "Submit"




- **Note:** The "Performed" section on the screen lets the user know how many racks have been harvested compared to the number of racks that have been planted & infected.
- The screen will return to the main menu

2. vHVT (Vaxxitek) Process

Below represents the example of the processes for the vHVT (Vaxxitek) virus.

- User will create the number of lot labels based on how many racks planted for vHVT process



- The screen will return to the main menu
- The labels will look like the following:






- After lot labels are printed, the user enters the Process control form:


```

C:\ Telnet 10.39.0.206
Merial MSCA AI Pro
1<Request lot label>
2<Process control
3<Edit control repo>
4<Exception process>
5<End of process >
6<Change Org >
7<Change Responsibi>
8<Logout >

```

- User associates all racks through the “Plant” process for the vHVT virus
- User will manually enter the cell start Lot number and plant number
- User presses  and chooses the lab they are currently working in
- User presses  and chooses the vHVT virus






- User presses  and the process will auto-populate with the “Plant”:

```


C:\ Telnet 10.39.0.206
Process AI
Lot Numbe>CEF01232B
Plant Num:JHV071621
Lab >TC1
Virus >vHVT (Uaxx
Process >1-Plant
Rack Code:
Ser Code :
<Continue>
<Submit>
<Cancel>

```

- User scans the Rack Code and Ser Code (lot) label
- User presses  to highlight “Continue” after each Rack and Ser Code is entered



- User presses  to highlight “Submit”



- All racks have now been associated with the Ser Code (lot) label under the “Plant” process.
- The screen will return to the main menu

- When user is ready to record the “Infect Seed” process for some of the racks, the user will go back into the Process control screen:

```

C:\ Telnet 10.39.0.206
Merial MSCA AI Pro
1<Request lot label>
2<Process control
3<Edit control repo>
4<Exception process>
5<End of process >
6<Change Org >
7<Change Responsibi>
8<Logout >


```

- User enter lot number, plant number, lab, and vHVT virus:

```

C:\ Telnet 10.39.0.206
Process AI
Lot Numbe>CEF01232B
Plant Num: JHV071621
Lab >IC1
Virus >JHVT <Uaxx
Process >
Rack Code:
Ser Code:
<Continue>
<Submit>
<Cancel>

```


- User presses  to choose the process
- Now that the racks are eligible for the next process, the user will now see the “Infect Seed” as a process option.
- **Note:** The “Plant” process is still an available option because more racks can be added if needed.


```


C:\ Telnet 10.39.0.206
Process
1<1-Plant >
2<2-Infect Seed>

```

- **Note:** User determines the Rack Codes to perform the “Infect Seed” process against. For this example, it will be for RACK1 and RACK2:



- User presses  and chooses “Submit” indicating that the “Infect Seed” process has been completed for these racks



- The screen will return to the main menu
- These racks have now been recorded in the system showing the “Infect Seed” process has been completed
- When user is ready to record the “Harvest Seed” process for the racks used in the “Infect Seed” process, the user will go back into the Process control screen:


```

C:\ Telnet 10.39.0.206
Merial MSCA AI Pro
1<Request lot label>
2<Process control
3<Edit control repo>
4<Exception process>
5<End of process >
6<Change Org >
7<Change Responsibi>
8<Logout >


```

- User enter lot number, plant number, lab, and vHVT virus:

```

C:\ Telnet 10.39.0.206
Process AI
Lot Numb>CEF01232B
Plant Num:UH0071621
Lab >IC1
Virus >vHVT (Uaxx
Process >
Rack Code:
Ser Code :
<Continue>
<Submit>
<Cancel>

```

- User presses  to choose the process
- Now that the racks are eligible for the next process, the user will now see the “Harvest Seed” as a process option.
- **Note:** The “Plant” and “Infect Seed” processes are still available options because more racks can be added if needed.

```

C:\ Telnet 10.39.0.206
Process
1<1-Plant >
2<2-Infect Seed >
3<3-Harvest Seed>

```

- User will enter the same racks for the “Harvest Seed” process. For this example, it will be for RACK1 and RACK2

Telnet 10.39.0.206

Process AI

Lot Numbe>CEF01232B

Plant Num:UHV071621

Lab >TC1

Virus >vHUT <Uaxx

Process >3-Harvest

Performed:

Rack Code:RACK1

Ser Code :35

<Continue>

<Submit>

<Cancel>

→

Telnet 10.39.0.206

Process AI

Lot Numbe>CEF01232B

Plant Num:UHV071621

Lab >TC1

Virus >vHUT <Uaxx

Process >3-Harvest

Performed:1 / 2

Rack Code:RACK2

Ser Code :36

<Continue>

<Submit>

<Cancel>

Telnet 10.39.0.206

Process AI

Lot Numbe>CEF01232B

Plant Num:UHV071621

Lab >TC1

Virus >vHUT <Uaxx

Process >3-Harvest

Performed:1 / 2

Rack Code:RACK2

Ser Code :36

<Continue>

<Submit>

<Cancel>

→

Telnet 10.39.0.206

Process AI

Lot Numbe>CEF01232B

Plant Num:UHV071621

Lab >TC1

Virus >vHUT <Uaxx

Process >3-Harvest

Performed:2 / 2

Rack Code:


Ser Code :


<Continue>

<Submit>

<Cancel>

- **Note:** Notice that the system remembers that there are only two racks eligible for the “Harvest Seed” process even though there are a total of four racks associated with this lot number in this example. (See the Performed section of the screen)

- User presses  and chooses “Submit” indicating that the “Harvest Seed” process has been completed for these racks




- The screen will return to the main menu
- These racks have now been recorded in the system showing the "Harvest Seed" process has been completed
- When the user is ready, the Exception process can be used to disassociate the racks used during the "Harvest Seed" process


G. Exception process for vHVT virus


To remove “Harvest Seed” processed racks:


- User will enter into the Exception process screen





- User will manually enter the lot number and plant number
- User will scan the Rack Code and Ser (lot) Code




- User presses  and a list of reason codes returns
- User chooses “Harvest Seed”



- User can choose to place a comment if desired, but is not mandatory



- User presses  to highlight "Save"
- The main menu returns
- User repeats the same exception process for the remaining "Harvest Seed" racks

```

C:\ Telnet 10.39.0.206
Exception Process AI
Lot Numbe>CEF01232B
Plant Num>UHV071621
Rack Code>RACK2
Ser Code >86
Reason >Harvest Se
Comment :
New Rack :
New Ser C:
<Save>
<Cancel>

```

- The main menu will return
- The “Harvest Seed” racks have now been disassociated.
- When user is ready to record the “Infect Serial” process for remaining racks, the user will go back into the Process control screen:

```

C:\ Telnet 10.39.0.206
Merial MSCA AI Pro
1<Request lot label>
2<Process control
3<Edit control repo>
4<Exception process>
5<End of process >
6<Change Org >
7<Change Responsibi>
8<Logout >

```

- User enter lot number, plant number, lab, and vHVT virus:

```

C:\ Telnet 10.39.0.206
Process AI
Lot Numbe>CEF01232B
Plant Num>UHV071621
Lab >IC1
Virus >vHVT (Vaxx
Process >
Rack Code:
Ser Code :
<Continue>
<Submit>
<Cancel>

```

- User presses  to choose the process


- Now that the racks are eligible for the next process, the user will now see the “Infect Serial” as a process option.
- Note:** The “Plant”, “Infect Seed”, “Harvest Seed” processes are still available options because more racks can be added if needed.

```

C:\ Telnet 10.39.0.206
Process
1<1-Plant      >
2<2-Infect Seed >
3<3-Harvest Seed >
4<4-Infect Serial>

```

- The remaining racks will be used from this point forward. For this example, it will be for RACK3 and RACK4



- User presses  to highlight “Submit”

```

C:\ Telnet 10.39.0.206
Process AI
Lot Numbe>CEF01232B
Plant Num:UHU071621
Lab      >TC1
Virus    >vHUT <Uaxx
Process  >4-Infect S
Performed:2 / 2
Rack Code:
Ser Code :
<Continue>
<Submit>
<Cancel>

```


- The main menu will return
- **Note:** Notice that there are only 2 out of 2 racks recognized for the “Infect Serial” process. This is because the previous racks used in the “Harvest Seed” process are no longer recognized because they have been disassociated through the Exception process.
- When user is ready to record the “Harvest Serial” process for remaining racks, the user will go back into the Process control screen:


```


C:\ Telnet 10.39.0.206
Merial MSCA AI Pro
1<Request lot label>
2<Process control
3<Edit control repo>
4<Exception process>
5<End of process   >
6<Change Org      >
7<Change Responsi>
8<Logout          >

```


- User enter lot number, plant number, lab, and vHVT virus:




- User presses  to choose the process
- Now that the racks are eligible for the next process, the user will now see the “Harvest Serial” as a process option.
- **Note:** The “Plant”, “Infect Seed”, “Harvest Seed”, and “Infect Serial” processes are still available options because more racks can be added if needed.



- The same racks used to complete the “Infect Serial” process will also be used in the “Harvest Serial” process. Again, for this example it will be for RACK3 and RACK4



- User presses  to highlight "Submit"




- The main menu will return


1. Exception Process:

The Exception Process is used when the technician has to remove the association between a rack, lot and plant number due to rack malfunction, contamination, Harvest Seed racks or other corrections prior to the end of the process. The technician can also provide the reason and any comments to assist in the explanation of the exception.

There are five reason codes that can be used during the exception process. However the "Harvest Seed" reason code will be available as an option only if the rack(s) entered have been recorded in the system as completed through the "Harvest Seed" process for the vHVT virus. Please see: Exception process for vHVT virus (removing "Harvest Seed" processed racks).



Versus



Important to remember:

Per the reason code chosen, the counts can be affected in different ways.

- "Clearing" will clear the prior process
- "Other" clears the racks in the count including the planted racks.
- "Unprocessed" will un-processes the chosen process.
- "Maintenance" will replace the rack with another. The data associated with the malfunctioning rack will be transferred to the replacement rack.

The Exception Process must be done before the Harvest process for all viruses except for the Harvest Seed process for the vHVT virus and the Unprocessed reason code; as the Harvest process is eligible for 'un-processing'.


A new Ser Code (lot) label does not have to be created if the original one can still be used.

H. Maintenance Reason Codes:

1. Maintenance Reason Code

A user will choose the 'Maintenance' reason code if the rack is not longer working properly and another rack will be used in its place. This reason code is used only if the user is going to replace the malfunctioning rack with another rack.


- **Note:** Before starting the exception process, determine if the lot label that is on the malfunctioning rack can be reused. If it cannot (i.e. it gets torn when trying to remove it from the rack or it won't re-stick securely to the new rack), then create a new lot label through the "Request lot label form". This new label will not need to be processed through process control form because the history will be transferred to the new label in the Exception process.
- User will choose the Exception process form



```

C:\ Telnet 10.39.0.206
Merial MSCA AI Pro
1<Request lot label>
2<Process control >
3<Edit control repo>
4<Exception process
5<End of process >
6<Change Org >
7<Change Responsibi>
8<Logout >

```

- User enter lot number and plant number
- User will scan the Rack Code and Ser Code
- User presses  to choose the reason




```

C:\ Telnet 10.39.0.206
Exception Process AI
Lot Numbe>CEF12141C
Plant Num>PLANT5
Rack Code>4321ABCD
Ser Code >28
Reason >
Comment :
New Rack :
New Ser C:
<Save>
<Cancel>

```


- User chooses "Maintenance"




- User can choose to place a comment if desired, but is not mandatory



- User scans the Rack Code which is on the replacement rack



- User presses  to highlight "Save"



- The main menu will return
- The following AI control report is generated through The “Edit control report” form and it is manually requested. (Please see Edit Control Report section for information on this report)
- Below is the reporting result of this exception

AI Control report

05-JUL-12

Page 1 / 1

Parameters Lot number : CEF12141C

Item : 24550000

Plant number : PLANT5

Combo Plant number :

Virus : HVT


Process : All

Virus :		HVT			
Process	Lab	Rack requested	Rack processed	Process date	Conclusion
1-Plant & Infect	TC1	5 (All)	5	03-JUL-12	Ok
		Rack code : 4321EFGH		Ser code : 29	
		Rack code : 5678ABCD		Ser code : 28	
2-Harvest	TC1	5 (All)	3	03-JUL-12	Count Does Not Match
		Rack code : 1234ABCD		Ser code : 25	
		Rack code : 1234EFGH		Ser code : 26	
		Rack code : 1234IJKL		Ser code : 27	

Exception list							
Rack code	Ser code	Reason	Comment	Date	User	Replaced by rack	Replaced by Ser
4321ABCD	28	Maintenance	Racks Stopped Turning	05-JUL-12 15:34:24	SIMSLETI	5678ABCD	28

2. Clearing Reason Code:

The “Clearing” reason code keeps the rack count in the last process completed; however, the rack will not be available to any remaining processes.

- User will choose the Exception process form
- User enter lot number and plant number
- User will scan the Rack Code and Ser Code
- User presses  to choose the reason

```

C:\ Telnet 10.39.0.206
Exception Process AI
Lot Numbe>CEF12141C
Plant Num>PLANT5
Rack Code>4321EFGH
Ser Code >29
Reason >
Comment :
New Rack :
New Ser C:
<Save>
<Cancel>


```

- User chooses “Clearing”

```

C:\ Telnet 10.39.0.206
Reason
1<Clearing>
2<Maintenance>
3<Other>
4<Unprocessed>

```

- User can choose to place a comment if desired, but is not mandatory
- User presses  to highlight “Save”

```

C:\ Telnet 10.39.0.206
Exception Process AI
Lot Numbe>CEF12141C
Plant Num>PLANT5
Rack Code>4321EFGH
Ser Code >29
Reason >Clearing
Comment :n/a
New Rack :
New Ser C:
<Save>
<Cancel>

```

- The main menu returns
- Below is the reporting result of this exception

AI Control report

05-JUL-12

Page 1 / 1

Parameters Lot number : CEF12141C

Item : 24550000

Plant number : PLANTS

Combo Plant number :

Virus : HVT


Process : All

Virus : HVT					
Process	Lab	Rack requested	Rack processed	Process date	Conclusion
1-Plant & Infect	TC1	5 (All)	5	03-JUL-12	Ok
		Rack code : 5678ABCD		Ser code : 28	
2-Harvest	TC1	4 (All)	3	03-JUL-12	Count Does Not Match
		Rack code : 1234ABCD		Ser code : 25	
		Rack code : 1234EFGH		Ser code : 26	
		Rack code : 1234IJKL		Ser code : 27	

Exception list							
Rack code	Ser code	Reason	Comment	Date	User	Replaced by rack	Replaced by Ser
4321ABCD	28	Maintenance	Racks Stopped Turning	05-JUL-12 15:34:24	SIMSLETI	5678ABCD	28
4321EFGH	29	Clearing (Plant & Infect)	n/a	05-JUL-12 16:05:08	SIMSLETI		

3. Other Reason Code:

“Other” reason code will remove a rack(s) from the count as if it never happened.


- User will choose the Exception process form
- User enter lot number and plant number
- User will scan the Rack Code and Ser Code
- User presses  to choose the reason

```

C:\ Telnet 10.39.0.206
Exception Process AI
Lot Numbe>CEF12141C
Plant Num>H071822
Rack Code>RACK26
Ser Code >57
Reason >-
Comment :
New Rack :
New Ser C:
<Save>
<Cancel>

```


- User chooses “Other”



- User can choose to place a comment if desired, but is not mandatory



- User presses  to highlight "Save"



- The main menu returns
- Below is the reporting result of this exception
- Report showing details before using the reason code "other":

AI Control report

18-JUL-12

Page 1 / 1

Parameters Lot number : CEF12141C

Item : 24550000

Plant number : H071822

Combo Plant number :

Virus : All

Process : All

Virus :		HVT			
Process	Lab	Rack requested	Rack processed	Process date	Conclusion
1-Plant & Infect	TC1	3 (All)	3	18-JUL-12	Ok
		Rack code : RACK26		Ser code : 57	
		Rack code : RACK27		Ser code : 58	
2-Harvest	TC1	3 (All)	1	18-JUL-12	Count Does Not Match
		Rack code : RACK25		Ser code : 56	

- Report showing details after using reason code "other":

AI Control report

18-JUL-12

Page 1 / 1

Parameters Lot number : CEF12141C

Item : 24550000

Plant number : H071822

Combo Plant number :

Virus : All


Process : All


Virus :		HVT			
Process	Lab	Rack requested	Rack processed	Process date	Conclusion
1-Plant & Infect	TC1, TC1	2 (All)	2	18-JUL-12	Ok
		Rack code : RACK27		Ser code : 58	
2-Harvest	TC1	2 (All)	1	18-JUL-12	Count Does Not Match
		Rack code : RACK25		Ser code : 56	

Exception list							
Rack code	Ser code	Reason	Comment	Date	User	Replaced by rack	Replaced by Ser
RACK26	57	Other		18-JUL-12 10:40:24	SIMSLETT		


4. Unprocessed Reason Code:

The “unprocessed” reason code is used when the user was to remove the rack from the “rack processed” column of the AI Process Report. It does not remove the rack from the “rack requested” column.


- User will choose the Exception process form
- User enter lot number, plant number,
- User will scan the Rack Code and Ser Code
- User presses  to choose the reason




- User chooses “Unprocessed”




- Choose the process for which the racks will be removed



- User can choose to place a comment if desired, but is not mandatory



- User presses  to highlight "Save"



- The main menu returns
- Below is the reporting result of this exception

Before using the "unprocessed" reason code:

AI Control report

18-JUL-12

Page 1 / 1

Parameters Lot number : CEF12141C

Item : 24550000

Plant number : H071822

Combo Plant number :

Virus : All

Process : All

Virus :		HVT			
Process	Lab	Rack requested	Rack processed	Process date	Conclusion
1-Plant & Infect	TC1	3 (All)	3	18-JUL-12	Ok
		Rack code : RACK26		Ser code : 57	
		Rack code : RACK27		Ser code : 58	
2-Harvest	TC1	3 (All)	1	18-JUL-12	Count Does Not Match
		Rack code : RACK25		Ser code : 56	

- After using the “unprocessed” reason code:

AI Control report

18-JUL-12

Page 1 / 1

Parameters Lot number : CEF12141C

Item : 24550000

Plant number : H071822

Combo Plant number :

Virus : All

Process : All

Virus :		HVT			
Process	Lab	Rack requested	Rack processed	Process date	Conclusion
1-Plant & Infect	TC1, TC1	2 (All)	2	18-JUL-12	Ok
		Rack code : RACK27		Ser code : 58	
2-Harvest	TC1	2 (All)	0	18-JUL-12	Count Does Not Match

Exception list							
Rack code	Ser code	Reason	Comment	Date	User	Replaced by rack	Replaced by Ser
RACK25	56	Unprocessed (Harvest Serial)	Contaminated	18-JUL-12 11:16:15	SIMSLETI		
RACK26	57	Other		18-JUL-12 11:11:04	SIMSLETI		


I. Edit Control Report

Edit control report → Call report (screen)


The Edit Control Report form provides the user with detail information of the events captured in the system associated with the cell start number and plant number. If you leave the virus and the process blank, you will get a separate report for each process completed. If you specify the virus and the plant, you will receive a report specifically for the virus and process requested.

Below are examples of how to create an AI Control Report:

- From the main menu, the user will choose “Edit control report”



- The user can enter only the lot and plant numbers, leaving other information blank, if desired



- The following report will return:

AI Control report

16-JUL-12

Page 1 / 1

Parameters Lot number : CEF01232B
 Plant number : VHV071621
 Combo Plant number :
 Virus : All
 Process : All

Item : 2455000B

Virus : vHVT (Vaxxitek)					
Process	Lab	Rack requested	Rack processed	Process date	Conclusion
1-Plant	TC1	4 (All)	4	16-JUL-12	Ok
2-Infect Seed	TC1	2	2	16-JUL-12	Ok
3-Harvest Seed	TC1	2	2	16-JUL-12	Ok
4-Infect Serial	TC1	2 (All)	2	16-JUL-12	Ok
5-Harvest Serial	TC1	2 (All)	2	16-JUL-12	Ok
		Rack code : RACK3		Ser code : 37	
		Rack code : RACK4		Ser code : 38	

Seed Racks Cleared					
Rack code	Ser code	Reason	Comment	Date	User
RACK1	35	Harvest Seed		16-JUL-12 14:24:47	SIMSLETI
RACK2	36	Harvest Seed		16-JUL-12 16:01:58	SIMSLETI

- The user can choose a specific virus

```

C:\ Telnet 10.39.0.206

Call report
Lot Numbe>CEF01232B
Plant Num>VHV071621
Combo Pla>
Virus    >vHVT <Vaxx
Process  >
<Submit>
<Cancel>

```

- The following report will return:

AI Control report

16-JUL-12
Page 1 / 1

Parameters Lot number : CEF01232B Item : 2455000B
 Plant number : VHV071621
 Combo Plant number :
 Virus : VHVT
 Process : All

Virus : VHVT (Vaxxitek)					
Process	Lab	Rack requested	Rack processed	Process date	Conclusion
1-Plant	TC1	4 (All)	4	16-JUL-12	Ok
2-Infect Seed	TC1	2	2	16-JUL-12	Ok
3-Harvest Seed	TC1	2	2	16-JUL-12	Ok
4-Infect Serial	TC1	2 (All)	2	16-JUL-12	Ok
5-Harvest Serial	TC1	2 (All)	2	16-JUL-12	Ok
		Rack code : RACK3		Ser code : 37	
		Rack code : RACK4		Ser code : 38	

Seed Racks Cleared					
Rack code	Ser code	Reason	Comment	Date	User
RACK1	35	Harvest Seed		16-JUL-12 14:24:47	SIMSLETI
RACK2	36	Harvest Seed		16-JUL-12 16:01:58	SIMSLETI

- The user can choose a specific virus **and** process

```

C:\ Telnet 10.39.0.206
Call report
Lot Numbe>CEF01232B
Plant Num>VHV071621
Combo Pla>
Virus      >VHVT <Vaxx
Process    >
<Submit>
<Cancel>

```



- User clicks  +  for the list of process values

```

C:\ Telnet 10.39.0.206
Process
1<1-Plant      >
2<2-Infect Seed >
3<3-Harvest Seed >
4<4-Infect Serial >
5<5-Harvest Serial >

```

```

C:\ Telnet 10.39.0.206
Call report
Lot Numbe>CEF01232B
Plant Num>VHV071621
Combo Pla>
Virus      >VHVT <Vaxx
Process    >3-Harvest
<Submit>
<Cancel>

```



- The following information is printed on the report:

AI Control report

16-JUL-12
Page 1 / 1

Parameters Lot number : CEF01232B Item : 2455000B
 Plant number : VHV071621
 Combo Plant number :
 Virus : VHVT
 Process : HARVEST_SEED

Virus : vHVT (Vaxxitek)					
Process	Lab	Rack requested	Rack processed	Process date	Conclusion
3-Harvest Seed	TC1	2	2	16-JUL-12	Ok

Seed Racks Cleared					
Rack code	Ser code	Reason	Comment	Date	User
RACK1	35	Harvest Seed		16-JUL-12 14:24:47	SIMSLETT
RACK2	36	Harvest Seed		16-JUL-12 16:01:58	SIMSLETT

For HVT / Rispens Combination processes, the “Combo Plant” feature can be used when printing a report. For this combination process, the lot number is the same but the Plant numbers will be different for each virus.

- User will enter the lot number and one of the plant numbers used
- User will enter the second plant number next to “Comb Pla(nt number)”

- User presses  to highlight Submit

```

C:\ Telnet 10.39.0.206
Call report
Lot Numbe>CEF01031B
Plant Num>MAINTTEST
Combo Pla>R071622
Virus >
Process >
<Submit>
<Cancel>
  
```

- The will receive two reports per virus
- The following **reports** will return

AI Control report

16-JUL-12

Page 1 / 2

Parameters Lot number : CEF01031B
 Plant number : MAINTTEST
 Combo Plant number : R071622
 Virus : All
 Process : All

Item : 2455000B

Virus : HVT					
Process	Lab	Rack requested	Rack processed	Process date	Conclusion
1-Plant & Infect	TC1	2 (All)	2	16-JUL-12	Ok
2-Harvest	TC1	2 (All)	2	16-JUL-12	Ok
		Rack code : RACK6		Ser code : 40	
		Rack code : RACK7		Ser code : 39	

Exception list							
Rack code	Ser code	Reason	Comment	Date	User	Replaced by rack	Replaced by Ser
RACK5	39	Maintenance		16-JUL-12 16:31:15	SIMSLETT	RACK7	39

AI Control report

16-JUL-12

Page 2 / 2


Parameters Lot number : CEF01031B
 Plant number : MAINTTEST
 Combo Plant number : R071622
 Virus : All
 Process : All

Item : 2455000B


Virus : Rispens					
Process	Lab	Rack requested	Rack processed	Process date	Conclusion
1-Plant & Infect	TC1	2 (All)	2	16-JUL-12	Ok
2-Harvest	TC1	2 (All)	2	16-JUL-12	Ok
		Rack code : RACK10		Ser code : 42	
		Rack code : RACK9		Ser code : 41	

J. Requesting AI Control Report from Oracle Form:

- Choose the **Merial Inventory US Select** responsibility
- Click View → Requests




- Click Submit a New Request...



The 'Find Requests' dialog box is shown with the following fields and options:

- Radio buttons for selection:
 - ☐ My Completed Requests
 - ☐ My Requests In Progress
 - ☐ All My Requests
 - ☒ Specific Requests
- Input fields for 'Specific Requests':
 - Request ID: []
 - Name: []
 - Date Submitted: []
 - Date Completed: []
 - Status: []
 - Phase: []
 - Requestor: []
- Checkbox: ☐ Include Request Set Stages in Query
- Order By: **Request ID** (dropdown menu)
- Select the Number of Days to View: **7** (text input)
- Buttons at the bottom: Submit a New Request..., Clear, Find


- Run a Single Request (defaulted)
- Click OK




The 'Submit a New Request' dialog box is shown with the following fields and options:

- Radio buttons for selection:
 - ☒ Single Request
 - ☐ Request Set
- Text descriptions:
 - Single Request: This allows you to submit an individual request.
 - Request Set: This allows you to submit a pre-defined set of requests.
- Buttons: OK, Cancel
- Select the Number of Days to View: **7** (text input)
- Buttons at the bottom: Submit a New Request..., Clear, Find

- Click the ... button for a list of reports
- Type %AI% and press Find
- Choose 'Merial Gainesville AI Control Report'



- Press OK
- Enter the Parameters (Note: the Lot Number and Plant Number must be entered)



- Click OK

Oracle Applications - BPTFOTE - Clone of BPTPRD as of 16-June-2012

File Edit View Folder Tools Window Help

Submit Request

Run this Request...

Copy...

Name: Merial Gainesville AI Control Report

Operating Unit:

Parameters: CEF01031B:R071622::RISPENS:PLANT_INFECT

Language: American English

Language Settings... Debug Options

At these Times...

Run the Job: As Soon as Possible

Schedule...

Upon Completion...

☒ Save all Output Files

Layout: Merial Gainesville AI Control Report

Options...

Notify:

Print to: Prod_TC3

Delivery Opts


Help (C) Submit Cancel


- User can see what printer the report will print to (Prod_TC3 = Lab TC3)
- Click Submit


K. End of Process


The “End of Process” allows the user to disassociate all racks used in a process so that the racks can be used for future virus processes. It is mandatory to print the AI Control Report prior to completing the “End of Process” in the system.


- From the main menu, the user chooses “End of process”






- User enters the lot number, plant number, lab, and virus
- Note: The virus will auto-populate when you press 



- **Note:** The user will need to manually print the report from the “Edit Control Report” form prior to entering the “End of Process” form. If the End of Process is done before the report is generated, there is no documentation available for reference afterwards.
- User presses  and will be required to answer if the “Report Pr(inted?)” by choosing “Yes” or “No”
- Choose “Yes”



- User presses  and Save is highlighted.
- Press  again



- The main menu returns
- the End of Process will disassociate all racks as expected

Note: If the user chooses “No” when asked if “Report Pr(inted?)”, then the system will not allow the End of Process to continue:

